

Product Datasheet

Recombinant Human Siglec-10 Protein, His Tag (orb1551243)

Description	Recombinant Human Siglec-10 Protein, His Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	Siglec10
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1244, Recombinant Human Mesothelin Protein is produced by mammalian expression system. The target protein is expressed with sequence (Met17-Thr546) of human Siglec-10 (Accession #Q96LC7) fused with a 6xHis Tag at the C-terminus.
Expiration Date	6 months from date of receipt.