

Product Datasheet

Recombinant Human CD3 epsilon & CD3 delta Heterodimer Protein, Fc His Tag (orb1551252)

Description	Recombinant Human CD3 epsilon & CD3 delta Heterodimer Protein, Fc His Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	CD3
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1235, Recombinant Human BAFFR Protein is produced by mammalian expression system. The target protein is expressed with sequence (Asp23-Asp126 (CD3E) & Phe22-Ala105 (CD3D)) of human CD3 epsilon & CD3 delta Heterodimer (Accession #P07766 & P04234) fused with an Fc tag at the C-terminus.
Expiration Date	6 months from date of receipt.