

Product Datasheet

Recombinant Human BTLA Protein, Avi & His Tag (orb1551258)

Description	Recombinant Human BTLA Protein, Avi & His Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	BTLA
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1229, Recombinant Human PD-L1 Protein is produced by mammalian expression system. The target protein is expressed with sequence (Met1-Ser150) of human BTLA (Accession #Q7Z6A9) fused with a 6xHis Tag, Avi tag at the C-terminus.
Expiration Date	6 months from date of receipt.