

## Product Datasheet

# Recombinant Human CD3 epsilon & CD3 gamma Heterodimer Protein, Fc His Tag (orb1551265)

<b>Description</b>	Recombinant Human CD3 epsilon & CD3 gamma Heterodimer Protein, Fc His Tag
<b>Species/Host</b>	Mammal
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 0.1 EU/ug by LAL method.
<b>Target</b>	CD3
<b>Note</b>	For research use only
<b>Purity</b>	> 95% by SDS-PAGE.
<b>Source</b>	KMP1222, Recombinant Human CD3 epsilon & CD3 delta Heterodimer Protein is produced by mammalian expression system. The target protein is expressed with sequence (Asp23-Asp126 (CD3E) & Gln23-Ser116 (CD3G)) of human CD3 epsilon & CD3 gamma Heterodimer (Accession #P07766 & P09693) fused with an Fc tag at the C-terminus.
<b>Expiration Date</b>	6 months from date of receipt.