

Product Datasheet

Biotinylated Recombinant Human NKG2A Protein, Avi & His Tag (orb1551266)

Description	Biotinylated Recombinant Human NKG2A Protein, Avi & His Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	NKG2A
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1221, Biotinylated Recombinant Human NKG2A Protein is produced by mammalian expression system. The target protein is expressed with sequence (Arg100-Leu233) of human NKG2A (Accession #P26715) fused with a 6xHis Tag, Avi tag at the N-terminus.
Expiration Date	6 months from date of receipt.