

## Product Datasheet

# Recombinant human FAP Protein, His Tag (orb1551281)

<b>Description</b>	Recombinant human FAP Protein, His Tag
<b>Species/Host</b>	Virus
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 0.1 EU/ug by LAL method.
<b>Target</b>	FAP
<b>Note</b>	For research use only
<b>Purity</b>	> 95% by SDS-PAGE.
<b>Source</b>	KMP1206, Recombinant human FAP Protein is produced by Baculovirus-Infected Sf9 Cells expression system. The target protein is expressed with sequence (Arg27-Asp760) of human FAP (Accession #NP_004451.2) fused with a 6xHis Tag at the C-terminus.
<b>Expiration Date</b>	6 months from date of receipt.