

Product Datasheet

Recombinant Human Osteoprotegerin/TNFRSF11B Protein, His Tag (orb1551284)

Description	Recombinant Human Osteoprotegerin/TNFRSF11B Protein, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	Osteoprotegerin
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1203, Recombinant Human Osteoprotegerin/TNFRSF11B Protein is produced by HEK293 expression system. The target protein is expressed with sequence (Glu22-Leu401) of human Osteoprotegerin/TNFRSF11B (Accession #NP_002537.3) fused with a 6xHis Tag at the C-terminus.
Expiration Date	6 months from date of receipt.