

Product Datasheet

Recombinant Human HER2/ErbB2 Protein, His Tag (orb1551287)

Description	Recombinant Human HER2/ErbB2 Protein, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	HER2
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1200, Recombinant Human HER2/ErbB2 Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Thr23-Thr652) of human HER2/ErbB2 (Accession #NP_004439.2) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.