

## Product Datasheet

# Recombinant Human PD-L1/B7-H1 Protein, Fc & His Tag (orb1551288)

<b>Description</b>	Recombinant Human PD-L1/B7-H1 Protein, Fc & His Tag
<b>Species/Host</b>	Human
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 0.1 EU/ug by LAL method.
<b>Target</b>	PDL1
<b>Note</b>	For research use only
<b>Purity</b>	> 95% by SDS-PAGE.
<b>Source</b>	KMP1199, Recombinant Human PD-L1/B7-H1 Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Phe19-Thr239) of human PD-L1/B7-H1 (Accession #NP_054862.1) fused with an Fc, 6×His Tag at the C-terminus.
<b>Expiration Date</b>	6 months from date of receipt.