

Product Datasheet

Recombinant Human CD30/TNFRSF8 Protein, His Tag (orb1551289)

Description	Recombinant Human CD30/TNFRSF8 Protein, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	CD30
Note	For research use only
Purity	> 97% by SDS-PAGE.
Source	KMP1198, Recombinant Human CD30/TNFRSF8 Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Phe19-Lys379) of human CD30/TNFRSF8 (Accession #NP_001234.2) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.