

Product Datasheet

Recombinant Human GFRA2/GFR α 2/GDNFRB Protein, His Tag (orb1551295)

Description	Recombinant Human GFRA2/GFR α 2/GDNFRB Protein, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	GFRA2
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1192, Recombinant Human GFRA2/GFR α 2/GDNFRB Protein is produced by HEK293 expression system. The target protein is expressed with sequence (Ser22-Ser441) of human GFRA2/GFR α 2/GDNFRB (Accession #NP_001486.4) fused with a 6 \times His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.