

## Product Datasheet

# Recombinant Human Contactin 2/CNTN2 Protein, His Tag (orb1551300)

<b>Description</b>	Recombinant Human Contactin 2/CNTN2 Protein, His Tag
<b>Species/Host</b>	Human
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 0.1 EU/ug by LAL method.
<b>Target</b>	Contactin 2
<b>Note</b>	For research use only
<b>Purity</b>	> 97% by SDS-PAGE.
<b>Source</b>	KMP1187, Recombinant Human Contactin 2/CNTN2 Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Ser31-Asn1012) of human Contactin 2/CNTN2 (Accession #NP_005067.1 ) fused with an 8×His Tag at the C-terminus.
<b>Expiration Date</b>	6 months from date of receipt.