

Product Datasheet

Recombinant Human SIRP alpha/CD172a Protein, Fc & His Tag (orb1551305)

Description	Recombinant Human SIRP alpha/CD172a Protein, Fc & His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	SIRP alpha
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1182, Recombinant Human SIRP alpha/CD172a Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Glu31-Arg370) of human SIRP alpha/CD172a (Accession #NP_542970.1) fused with an Fc, 6xHis Tag at the C-terminus.
Expiration Date	6 months from date of receipt.