

Product Datasheet

Recombinant Human IL-21R/Interleukin-21 Receptor Protein, His Tag (orb1551309)

Description	Recombinant Human IL-21R/Interleukin-21 Receptor Protein, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	IL21R
Note	For research use only
Purity	> 90% by SDS-PAGE.
Source	KMP1178, Recombinant Human IL-21R/Interleukin-21 Receptor Protein is produced by HEK293 expression system. The target protein is expressed with sequence (Cys20-Pro236) of human IL-21R/Interleukin-21 Receptor (Accession #NP_068570.1) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.