

Product Datasheet

Recombinant Human PDGFRa/CD140a Protein, Fc & His Tag (orb1551316)

Description	Recombinant Human PDGFRa/CD140a Protein, Fc & His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	PDGFRa
Note	For research use only
Purity	> 97% by SDS-PAGE.
Source	KMP1171, Recombinant Human PDGFRa/CD140a Protein is produced by HEK293 expression system. The target protein is expressed with sequence (Gln24-Glu524) of human PDGFRa/CD140a (Accession #NP_006197.1) fused with an Fc, 6xHis Tag at the C-terminus.
Expiration Date	6 months from date of receipt.