

Product Datasheet

Recombinant Human LDLR/LDL Receptor Protein, His Tag (orb1551317)

Description	Recombinant Human LDLR/LDL Receptor Protein, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	LDLR
Note	For research use only
Purity	> 90% by SDS-PAGE.
Source	KMP1170, Recombinant Human LDLR/LDL Receptor Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Ala22-Arg788) of human LDLR/LDL Receptor (Accession #NP_000518.1) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.