

Product Datasheet

Recombinant Human KLK11/Kallikrein-11 Protein, His Tag (orb1551335)

Description	Recombinant Human KLK11/Kallikrein-11 Protein, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	KLK11
Note	For research use only
Purity	> 92% by SDS-PAGE.
Source	KMP1152, Recombinant Human KLK11/Kallikrein-11 Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Met1-Asn250) of human KLK11/Kallikrein-11 (Accession #NP_006844.1) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.