

Product Datasheet

Recombinant Human Ephrin-A4/EFNA4 Protein, Fc & His Tag (orb1551342)

Description	Recombinant Human Ephrin-A4/EFNA4 Protein, Fc & His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	EphrinA4
Note	For research use only
Purity	> 80% by SDS-PAGE.
Source	KMP1145, Recombinant Human Ephrin-A4/EFNA4 Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Leu26-Gly171) of human Ephrin-A4/EFNA4 (Accession # NP_005218.1) fused with an Fc, 6xHis Tag at the C-terminus.
Expiration Date	6 months from date of receipt.