

Product Datasheet

Recombinant Human Mucin-1/MUC-1 Protein, Fc & His Tag (orb1551344)

Description	Recombinant Human Mucin-1/MUC-1 Protein, Fc & His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	Mucin1
Note	For research use only
Purity	> 97% by SDS-PAGE.
Source	KMP1143, Recombinant Human Mucin-1/MUC-1 Protein is produced by HEK293 expression system. The target protein is expressed with sequence (Met1-Gly167) of human Mucin-1/MUC-1 (Accession #NP_001018016.1) fused with an Fc, 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.