

## Product Datasheet

# Recombinant Human PRLR/Prolactin Receptor Protein, His Tag (orb1551348)

|                        |  |
|------------------------|--|
| <b>Description</b>     | Recombinant Human PRLR/Prolactin Receptor Protein, His Tag   |
| <b>Species/Host</b>    | Human  |
| <b>Reactivity</b>      | Human  |
| <b>Endotoxins</b>      | < 0.1 EU/ug by LAL method.   |
| <b>Target</b>          | PRLR   |
| <b>Note</b>            | For research use only  |
| <b>Purity</b>          | > 97% by SDS-PAGE.   |
| <b>Source</b>          | KMP1139, Recombinant Human PRLR/Prolactin Receptor Protein is produced by HEK293 expression system. The target protein is expressed with sequence (Gln25-Asp234) of human PRLR/Prolactin Receptor (Accession #NP_000940.1) fused with a 6xHis Tag at the C-terminus. |
| <b>Expiration Date</b> | 6 months from date of receipt.   |