

## Product Datasheet

# Recombinant Human GFRA1/GFR $\alpha$ 1/GDNFRA Protein, Fc & His Tag (orb1551356)

<b>Description</b>	Recombinant Human GFRA1/GFR $\alpha$ 1/GDNFRA Protein, Fc & His Tag
<b>Species/Host</b>	Human
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 1.0 EU/ug by LAL method.
<b>Target</b>	GFRA1
<b>Note</b>	For research use only
<b>Purity</b>	> 80% by SDS-PAGE.
<b>Source</b>	KMP1131, Recombinant Human GFRA1/GFR $\alpha$ 1/GDNFRA Protein is produced by HEK293 expression system. The target protein is expressed with sequence (Ala19-Lys429) of human GFRA1/GFR $\alpha$ 1/GDNFRA (Accession #NP_665736) fused with an Fc, 6xHis Tag at the C-terminus.
<b>Expiration Date</b>	6 months from date of receipt.