

Product Datasheet

Recombinant Human Tie2/CD202b/TEK Protein, Fc & His Tag (orb1551360)

Description	Recombinant Human Tie2/CD202b/TEK Protein, Fc & His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	Tie2
Note	For research use only
Purity	> 97% by SDS-PAGE.
Source	KMP1127, Recombinant Human Tie2/CD202b/TEK Protein is produced by HEK293 expression system. The target protein is expressed with sequence (Ala23-Lys745) of human Tie2/CD202b/TEK (Accession #NP_000450.2) fused with an Fc, 6xHis Tag at the C-terminus.
Expiration Date	6 months from date of receipt.