

Product Datasheet

Recombinant Human ALCAM/CD166 Protein, Fc & His Tag (orb1551362)

Description	Recombinant Human ALCAM/CD166 Protein, Fc & His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	ALCAM
Note	For research use only
Purity	> 97% by SDS-PAGE.
Source	KMP1125, Recombinant Human ALCAM/CD166 Protein is produced by HEK293 expression system. The target protein is expressed with sequence (Met1-Ala526) of human ALCAM/CD166 (Accession #NP_001618.2) fused with an Fc, 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.