

Product Datasheet

Recombinant Human SFTPD/SP-D/SFTP4 Protein, His Tag (orb1551363)

Description	Recombinant Human SFTPD/SP-D/SFTP4 Protein, His Tag
Reactivity	Human
Endotoxins	> 97% by SDS-PAGE.
Target	SFTPD
Note	For research use only
Purity	> 97% by SDS-PAGE.
Source	Recombinant Human SFTPD/SP-D/SFTP4 Protein is produced by HEK293 expression system. The target protein is expressed with sequence (Ala21-Phe375 (Glu22Gly)) of human SFTPD/SP-D/SFTP4 (Accession #NP_003010.4) fused with a 6His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)