

Product Datasheet

Recombinant Human Cystatin C/CST3 Protein, His Tag (orb1551368)

Description	Recombinant Human Cystatin C/CST3 Protein, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	Cystatin C
Note	For research use only
Purity	> 97% by SDS-PAGE.
Source	KMP1119, Recombinant Human Cystatin C/CST3 Protein is produced by HEK293 expression system. The target protein is expressed with sequence (Met1-Ala146) of human Cystatin C/CST3 (Accession #NP_000090.1) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.