

Product Datasheet

Recombinant Human EPO Receptor/EPOR Protein, Fc & His Tag (orb1551369)

Description	Recombinant Human EPO Receptor/EPOR Protein, Fc & His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	EPO
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1118, Recombinant Human EPO Receptor/EPOR Protein is produced by HEK293 expression system. The target protein is expressed with sequence (Ala25-Pro250) of human EPO Receptor/EPOR (Accession #NP_000112.1) fused with an Fc, 6xHis Tag at the C-terminus.
Expiration Date	6 months from date of receipt.