

Product Datasheet

Recombinant Human Cystatin SN/CST1 Protein, His Tag (orb1551374)

Description	Recombinant Human Cystatin SN/CST1 Protein, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	Cystatin SN
Note	For research use only
Purity	> 97% by SDS-PAGE.
Source	KMP1113, Recombinant Human Cystatin SN/CST1 Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Trp21-Ser141) of human Cystatin SN/CST1 (Accession #NP_001889.2) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.