

Product Datasheet

Recombinant Human VEGF R2/KDR Protein, His Tag (orb1551384)

Description	Recombinant Human VEGF R2/KDR Protein, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	VEGF R2
Note	For research use only
Purity	> 97% by SDS-PAGE.
Source	KMP1103, Recombinant Human VEGF R2/KDR Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Ala20-Glu764) of human VEGF R2/KDR (Accession #NP_002244.1) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.