

Product Datasheet

Recombinant Human MANF Protein (orb1551391)

Description	Recombinant Human MANF Protein
Species/Host	E. coli
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	MANF
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1096, Recombinant Human MANF Protein is produced by E. coli expression system. The target protein is expressed with sequence (Leu25-Leu182) of human MANF (Accession #NP_006001.4).
Expiration Date	6 months from date of receipt.