

Product Datasheet

Recombinant Human FGF-12 Protein, His Tag (orb1551397)

Description	Recombinant Human FGF-12 Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	< 1.0 EU/ug by LAL method.
Target	FGF12
Note	For research use only
Purity	> 97% by SDS-PAGE.
Source	KMP1090, Recombinant Human FGF-12 Protein is produced by E. coli expression system. The target protein is expressed with sequence (Met1-Thr181) of human FGF-12 (Accession #NP_004104.3) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.