

Product Datasheet

Recombinant Human ErbB3/Her3 Protein, His Tag (orb1551401)

Description	Recombinant Human ErbB3/Her3 Protein, His Tag
Species/Host	Human
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	ErbB3
Note	For research use only
Purity	> 97% by SDS-PAGE.
Source	KMP1086, Recombinant Human ErbB3/Her3 Protein is produced by HEK293 expression system. The target protein is expressed with sequence (Ser20-Thr643) of human ErbB3/Her3 (Accession #NP_001973.2) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.