

Product Datasheet

Recombinant Human MCP-3/CCL7 Protein (orb1551406)

Description	Recombinant Human MCP-3/CCL7 Protein
Species/Host	E. coli
Reactivity	Human
Endotoxins	< 1.0 EU/ug by LAL method.
Target	MCP3
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1081, Recombinant Human MCP-3/CCL7 Protein is produced by E. coli expression system. The target protein is expressed with sequence (Gln24-Leu99) of human MCP-3/CCL7 (Accession #NP_006264.2).
Expiration Date	6 months from date of receipt.