

Product Datasheet

Recombinant Human BCL2 Protein (orb1551410)

Description	Recombinant Human BCL2 Protein
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact us for more information.
Target	BCL2
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1077, Recombinant Human BCL2 Protein is produced by E. coli expression system. The target protein is expressed with sequence (Ala2-Asp211) of human BCL2 (Accession #NP_000624.2).
Expiration Date	6 months from date of receipt.