

## Product Datasheet

# Recombinant Human Galectin-3 Protein (orb1551411)

<b>Description</b>	Recombinant Human Galectin-3 Protein
<b>Species/Host</b>	E. coli
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 0.1 EU/ug by LAL method.
<b>Target</b>	Galectin3
<b>Note</b>	For research use only
<b>Purity</b>	> 97% by SDS-PAGE.
<b>Source</b>	KMP1076, Recombinant Human Galectin-3 Protein is produced by E. coli expression system. The target protein is expressed with sequence (Ala2-Ile250) of human Galectin-3 (Accession #NP_002297.2).
<b>Expiration Date</b>	6 months from date of receipt.