

## Product Datasheet

### Recombinant Human cAspAT Protein (orb1551412)

<b>Description</b>	Recombinant Human cAspAT Protein
<b>Species/Host</b>	E. coli
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 0.1 EU/ug by LAL method.
<b>Target</b>	cAspAT
<b>Note</b>	For research use only
<b>Purity</b>	> 95% by SDS-PAGE.
<b>Source</b>	KMP1075, Recombinant Human cAspAT Protein is produced by E. coli expression system. The target protein is expressed with sequence (Ala2-Gln413) of human cAspAT (Accession #NP_002070.1).
<b>Expiration Date</b>	6 months from date of receipt.