

## Product Datasheet

# Recombinant Human IL-1RN Protein, His Tag (orb1551413)

<b>Description</b>	Recombinant Human IL-1RN Protein, His Tag
<b>Species/Host</b>	E. coli
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 0.1 EU/ug by LAL method.
<b>Target</b>	IL1RN
<b>Note</b>	For research use only
<b>Purity</b>	> 97% by SDS-PAGE.
<b>Source</b>	KMP1074, Recombinant Human IL-1RN Protein is produced by E. coli expression system. The target protein is expressed with sequence (Arg26-Glu177) of human IL-1RN (Accession #NP_776214.1) fused with an initial Met at the N-terminus and a 6×His Tag at the C-terminus.
<b>Expiration Date</b>	6 months from date of receipt.