

Product Datasheet

Recombinant Human IL-33 Protein (orb1551414)

Description	Recombinant Human IL-33 Protein
Species/Host	E. coli
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	IL33
Note	For research use only
Purity	> 97% by SDS-PAGE.
Source	KMP1073, Recombinant Human IL-33 Protein is produced by E. coli expression system. The target protein is expressed with sequence (Ser112-Thr270) of human IL-33 (Accession #NP_001300973.1).
Expiration Date	6 months from date of receipt.