

Product Datasheet

Recombinant Human 4-1BB Ligand/TNFSF9 Protein, His Tag (orb1551415)

Description	Recombinant Human 4-1BB Ligand/TNFSF9 Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	41BB
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1072, Recombinant Human 4-1BB Ligand/TNFSF9 Protein is produced by E. coli expression system. The target protein is expressed with sequence (Arg71-Glu254) of human 4-1BB Ligand/TNFSF9 (Accession #NP_003802.1) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.