

Product Datasheet

Recombinant Human Leptin Protein (orb1551418)

Description	Recombinant Human Leptin Protein
Species/Host	E. coli
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	Leptin
Note	For research use only
Purity	> 92% by SDS-PAGE.
Source	KMP1069, Recombinant Human Leptin Protein is produced by E. coli expression system. The target protein is expressed with sequence (Val22-Cys167) of human Leptin (Accession #NP_000221.1).
Expiration Date	6 months from date of receipt.