

Product Datasheet

Recombinant Human IL-16 Protein, His Tag (orb1551419)

Description	Recombinant Human IL-16 Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	IL16
Note	For research use only
Purity	> 97% by SDS-PAGE.
Source	KMP1068, Recombinant Human IL-16 Protein is produced by E. coli expression system. The target protein is expressed with sequence (Pro2- Ser130) of human IL-16 (Accession #NP_757366.2) fused with a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.