

Product Datasheet

Recombinant Human FKBP12 Protein (orb1551420)

Description	Recombinant Human FKBP12 Protein
Species/Host	E. coli
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	FKBP12
Note	For research use only
Purity	> 90% by SDS-PAGE.
Source	KMP1067, Recombinant Human FKBP12 Protein is produced by E. coli expression system. The target protein is expressed with sequence (Met1-Glu108) of human FKBP12 (Accession #NP_000792.1).
Expiration Date	6 months from date of receipt.