

Product Datasheet

Recombinant Human IL-11 Protein (orb1551422)

Description	Recombinant Human IL-11 Protein
Species/Host	E. coli
Reactivity	Human
Endotoxins	< 1.0 EU/ug by LAL method.
Target	IL11
Note	For research use only
Purity	> 97% by SDS-PAGE.
Source	KMP1065, Recombinant Human IL-11 Protein is produced by E. coli expression system. The target protein is expressed with sequence (Pro22-Leu199) of human IL-11 (Accession #NP_000632.1) fused with an initial Met at the N-terminus.
Expiration Date	6 months from date of receipt.