

## Product Datasheet

### Recombinant Human PIN1 Protein (orb1551423)

<b>Description</b>	Recombinant Human PIN1 Protein
<b>Species/Host</b>	E. coli
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 0.1 EU/ug by LAL method.
<b>Target</b>	PIN1
<b>Note</b>	For research use only
<b>Purity</b>	> 95% by SDS-PAGE.
<b>Source</b>	KMP1064, Recombinant Human PIN1 Protein is produced by E. coli expression system. The target protein is expressed with sequence (Met1-Glu163) of human PIN1 (Accession #NP_006212.1).
<b>Expiration Date</b>	6 months from date of receipt.