

Product Datasheet

Recombinant Human FGFR-2 Protein (orb1551425)

Description	Recombinant Human FGFR-2 Protein
Species/Host	Virus
Reactivity	Human
Endotoxins	< 1.0 EU/ug by LAL method.
Target	FGFR2
Note	For research use only
Purity	> 97% by SDS-PAGE.
Source	KMP1062, Recombinant Human FGFR-2 Protein is produced by insect Cell-baculovirus expression system. The target protein is expressed with sequence (Pro458-Gln778) of human FGFR-2 (Accession #NP_000132.3).
Expiration Date	6 months from date of receipt.