

## Product Datasheet

### Recombinant Human BCL2 Protein (orb1551426)

<b>Description</b>	Recombinant Human BCL2 Protein
<b>Species/Host</b>	E. coli
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 1.0 EU/ug by LAL method.
<b>Target</b>	BCL2
<b>Note</b>	For research use only
<b>Purity</b>	> 90% by SDS-PAGE.
<b>Source</b>	KMP1061, Recombinant Human BCL2 Protein is produced by E. coli expression system. The target protein is expressed with sequence (Ala2-Asp211) of human BCL2 (Accession #NP_000624.2) fused with no tags.
<b>Expiration Date</b>	6 months from date of receipt.