

Product Datasheet

Recombinant Human Sonic hedgehog N-product Protein (orb1551427)

Description	Recombinant Human Sonic hedgehog N-product Protein
Species/Host	E. coli
Reactivity	Human
Endotoxins	< 1.0 EU/ug by LAL method.
Target	Sonic hedgehog Nproduct
Note	For research use only
Purity	> 97% by SDS-PAGE.
Source	KMP1060, Recombinant Human Sonic hedgehog N-product Protein is produced by E. coli expression system. The target protein is expressed with sequence (Cys24-Gly197 (Cys24Ile-Val-Ile)) of human Sonic hedgehog N-product (Accession #NP_000184.1).
Expiration Date	6 months from date of receipt.