

Product Datasheet

Recombinant Human ISG15 Protein (orb1551428)

Description	Recombinant Human ISG15 Protein
Species/Host	E. coli
Reactivity	Human
Endotoxins	< 1.0 EU/ug by LAL method.
Target	ISG15
Note	For research use only
Purity	> 97% by SDS-PAGE.
Source	KMP1059, Recombinant Human ISG15 Protein is produced by E. coli expression system. The target protein is expressed with sequence (Gly2-Gly157) of human ISG15 (Accession #NP_005092.1).
Expiration Date	6 months from date of receipt.