

Product Datasheet

Recombinant Human NRG1-alpha/HRG1-alpha Protein (orb1551430)

Description	Recombinant Human NRG1-alpha/HRG1-alpha Protein
Species/Host	E. coli
Reactivity	Human
Endotoxins	< 1.0 EU/ug by LAL method.
Target	NRG1alpha
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP1057, Recombinant Human NRG1-alpha/HRG1-alpha Protein is produced by E. coli expression system. The target protein is expressed with sequence (Ser177-Lys241) of human NRG1-alpha/HRG1-alpha (Accession #NP_039258.1).
Expiration Date	6 months from date of receipt.