

Product Datasheet

Recombinant Human LRPAP Protein, His Tag (orb1551432)

Description	Recombinant Human LRPAP Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	< 1.0 EU/ug by LAL method.
Target	LRPAP
Note	For research use only
Purity	> 87% by SDS-PAGE.
Source	KMP1055, Recombinant Human LRPAP Protein is produced by E. coli expression system. The target protein is expressed with sequence (Tyr35-Leu357) of human LRPAP (Accession #NP_002328.1) fused with an initial Met at the N-terminus and a 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.