

## Product Datasheet

# Recombinant Human NEDD8 Protein, His Tag (orb1551433)

<b>Description</b>	Recombinant Human NEDD8 Protein, His Tag
<b>Species/Host</b>	E. coli
<b>Reactivity</b>	Human
<b>Endotoxins</b>	Please contact us for more information.
<b>Target</b>	NEDD8
<b>Note</b>	For research use only
<b>Purity</b>	> 97% by SDS-PAGE.
<b>Source</b>	KMP1054, Recombinant Human NEDD8 Protein is produced by E. coli expression system. The target protein is expressed with sequence (Met1-Gly76) of human NEDD8 (Accession #NP_006147.1) fused with a 6×His Tag at the C-terminus.
<b>Expiration Date</b>	6 months from date of receipt.